

The Role of Reaction Calorimetry in the Development of Safe, Efficient and Reliable Chemical Processes



Agenda

1	What is Reaction Calorimetry?
2	Reasons to Perform Reaction Calorimetry
3	Technical Aspects
4	Our Reaction Calorimetry Market Offering
5	Appendix



Definitions& Concepts

In general, calorimetry is the measurement of heat.

Heat is thermal energy, which derives from the vibrations and movement of molecules Heat is usually released when bonds are formed between atoms

The rate of heat release is proportional to the overall rate of the reaction

In the first instance, releasing heat increases the temperature of the reaction mixture

Calorimetry requires accurate control and measurement of temperature

METTLER TOLEDO AutoChem provides solutions for Reaction Calorimetry

Reaction calorimetry is a non-destructive, real-time technique that provides chemical process data. It is widely used by process and development chemists as well as dedicated process safety engineers. A reaction calorimeter accurately measures the enthalpy changes of chemical reactions ("reaction enthalpy")



When does a reaction start?

- When does a reaction stop?
- What factors influence the reaction rate?

How much energy is released, and when?

- Calorimetry can also help to derive:
- Reaction mechanisms, pathways
- Reaction Kinetics
- Critical Process Parameters, CPP (dosing, mixing ...)
- Scale up parameters (KJ/kg, cooling potential ...)
- Safety Parameters (ATad, MTSR ...)





Reaction Enthalpies vary widely between reactions. It is obvious that the amount of energy released is directly linked to the potential damage in case of an accident.

Typical Example of	Energy	
Reaction Type	H kJ.mol ⁻¹	
Neutralization (HCI)	-55	
Hydrolysis (Acetic Anhydride)	-57	Exothermic
Polymerization (Styrene)	-60	
Epoxidation	-100	
Neutralization (H2SO4)	-105	
Amination	-120	
Nitration	-130	
Sulphonation	-150	
Hydrogenation (Alkene)	-200	Strongly
Hydrogenation (Nitro)	-560	Exothermic
	KODUTAL CANERA OD VOL	


ROCAELI SANAYI ODASI PR SES EMNIYETI SEMPOZYUMU

Why is Reaction Calorimetry Critical?

Reaction Calorimetry provides the basic information to create safe and scalable processes


Provides information about progress reaction

- Induction time
- Start and End of reaction
- Maximum heat flow
- Reaction enthalpy
- Heat flow in function of dosing



Provides information about potential safety issues


- Accumulation of reactants
- Non-scalable conditions





Where and When is Reaction Calorimetry Applied?


The Three Pillars of Application





MT calorimeters are optimized primarily for liquid-phase reactions


The focus on process safety, scale up and process understanding (see _____)






Reasons to Perform Calorimetry in R&D

1. Reaction / Process Understanding – Points of Interest



KOCAELİ SANAYİ ODASI PR SES

Reaction and Process Understanding



- Check behaviour of process at small scale (100ml or less)
- Assess scalability (heat release)
- Identify critical process steps
- Rapid screening of process parameters
- Parallel reactions




Reasons to Perform Calorimetry: Scaleup

2. Process Scale-Up – Points of Interest

Energy balance Temperature Profile Cooling Requirement Mixing / Mass Transfer Dosing profile Yield / Quality Optimization Fit to Plant





"Got a few problems going from lab scale up to full-scale commercial."

ROCAELI SANAYI ODASI

Process Scale– Up





- Assessment of Enthalpy, Heat Capacity, Heat Transfer & Mass Transfer under Process-Like Conditions (up to 1 litre scale)
- Develop dosing and temperature profiles
- Optimize critical process steps



Reasons to Perform Calorimetry in Safety

3. Process Safety – Points of Interest

Ensure safety of process Understand limits Evaluate risks Avoid disasters Protect the environment Save lives




Magnablend, Waxahachie TX, October 2011



KOCAELİ SANAYİ ODASI PRUSESES EMNIYETI SEMPOZYUMU


Safety Concept Assessment



PR SES

Why Perform Risk Assessment?

The 6 questions that help to develop the runaway scenario and provide guidance for the determination of data required for the risk assessment:





Can the process temperature be controlled by the cooling system?



Question 2

What temperature can be attained after runaway of the desired reaction?



Question 3

What temperature can be attained after runaway of the secondary reaction?



Question 4

moment does


failure have the

consequences?

the cooling

At which

worst



Question 5

How fast is the runaway of the desired reaction?




Question 6

How fast is the runaway of the decomposition starting at MTSR?

KOCAELİ SANAYİ ODASI PR SES EMNIYETI SEMPOZYUMU

Process Safety



- Assessment of Enthalpy, Heat Capacity, Heat Transfer, Mass Transfer under Process-Like Conditions
- Safely explore wide range of conditions and process upsets
- Quantify the process risk

KOCAELİ SANAYİ ODASI PR SES

Bad Calorimetry Data = Lack of understanding of thermochemistry Potential risk for accidents

For the User

- Difficulty interpreting the results equates to a waste of time
- Inaccurate results
- Risk of personnel injuries

For the Management

- No trust in data Plant cannot be run in an optimal way
- Need to apply large "safety buffers" resulting in a less efficient process
- Risk of accidents

For the Company

- Environmental consequences
- Monetary and reputation loss



What Good and Reliable Calorimetry Data Means

Accurate and Precise results will allow you to design a better and safer process at lower cost (when you need it most!)

- An integrated thermostat will
- have better and faster control
- ensure that system can be run isothermally
- Provide heat capacity (Cp) data and good results under non-isothermal conditions
- be safer in case of emergency

"Simple experiments" will work with most systems, however ...

Calorimetry with external circulator/cryostat will be more difficult if the reaction

- is highly exothermic
- is non-isothermal
- has large changes in Heat Transfer properties (viscosity, precipitation etc.)

For challenging experiments, the accuracy of the RC1mx is much better

The RC1mx will get good results for reactions under difficult but real conditions

- close to boiling temperature
- under reflux



Our Reaction Calorimetry Market Offering

Mettler-Toledo AutoChem provides you with:

Global Support and Service with local contacts Over 30 years of Applications Know-How and Experience Over 1000 satisfied users across the chemical industry Web support, User Meetings, International Conferences AutoChem Community Website A comprehensive range of options and accessories Continuous research and development of new technologies Unique RTCal Real Time Calorimetry – save time and money Modern, reliable hardware and software – maximum utilization rates Proven ROI – case studies available





